

Differential Scanning Calorimetry (DSC)



Q-series
Discovery
DSC2500
DSC250
DSC25



©2018 TA Instruments

Calibration & Verification



The DSC Heat Flow Rate Equation

- A DSC measures the difference in Heat Flow Rate between a sample and inert reference as a function of time and temperature.

$$\frac{dH}{dt} = C_p \frac{dT}{dt} + f(T, t)$$

- A DSC is calibrated for the heat flow enthalpy and temperature. Baseline calibrations are performed per manufacturers recommendations.

Calibration of Specific Instrument Models

Tzero

Measurement of
Rs and Cs

DSC2500

DSC250

1st Generation
Discovery

Q2000

Q200

Baseline calibration



Measurement of
slope and offset

DSC 25

Q20

Cell constant and temperature

All DSCs



General Calibration and Verification Guidelines

- Calibration
 - Use Calibration Mode
 - Calibrate upon installation
 - Re-calibrate if does not pass verification or if instrument setup is modified (see previous slide)
- Verification
 - Determine how often to verify data
 - Run a reference material as a sample (in standard mode)
 - Compare results vs literature values
 - If results are within your tolerance – system checks out and does not need re-calibration
 - If results are out of tolerance, then re-calibrate

Requirements Prior to Calibration

- DSC cell must be free of contaminants
- An inert purge gas, such as nitrogen, where the flow rate is controlled to 10-50 ml/min +/- 5 ml/min
- A balance to weigh specimens and containers to at least +/- 0.1 mg. The balance should have a capacity greater than 20 mg.
- High purity reference materials (>99.99%) for calibration

Instrument Setup Factors Affecting Calibration

- Purge Gas Type
 - Re-calibrate baseline/Tzero, temperature and cell constant
 - Thermal conductivity of helium \neq Thermal conductivity of nitrogen/air/oxygen \neq Thermal conductivity of argon
- Cooling Accessories
 - Re-calibrate baseline/Tzero, temperature and cell constant
 - The position of the cooling head around the cell will affect the calibration of the instrument. Uninstallation and reinstallation of a cooling accessory or changing the cooling accessory warrants a complete re-calibration
- Pan Selection
 - Re-calibrate temperature and cell constant
 - It will not impact the baseline/Tzero calibration

ASTM E 967 - Standard Test Method for Temperature Calibration of DSC's

- For transition temperature, calibration is required with known reference standards.
 - Pure metals (indium, lead, etc.) typically used
 - Extrapolated onset is used as melting temperature
 - Sample is fully melted at the peak
- This test method consists of heating the calibration materials at a controlled rate in a controlled atmosphere through a region of known thermal transition.
- The heat flow into the calibration material or the difference of temperature between the calibration material and a reference is monitored and continuously recorded.
- A transition is marked by the absorption of energy by the specimen resulting in a corresponding endothermic peak in the heating curve.

ASTM E 967 - Temperature Calibration of DSC's



Indium (156.6°C)
Can be re-used, can pre-melt



Tin (231.9°C)
Has multiple crystalline forms,
use once – no pre-melt




Zinc (419.5°C)
Can easily oxidize and alloy with
container, use once – no pre-melt

ASTM E 968 - Standard Practice for Heat Flow (Enthalpy) Calibration of DSC's

- For enthalpy of a transition, calibration is required with known reference standards.
- This standard consists of calibrating the heat flow response of a DSC by recording the melting endotherm of a high-purity standard material as a function of time.
- The peak is then integrated (over time) to yield an area measurement proportional to the enthalpy of melting of the standard material.
- Indium is a commonly used reference standard.

Temperature and Cell Constant Calibration


- Prepare a 3-5 mg sample of indium and “pre-melt” prior to first use



- Verify at least once a month
- Typical values for cell constant:
 - 0.9 to 1.2 (in N₂)
 - Helium will typically give higher values for the cell constant

Calorimetric and Temperature Calibration (Un-calibrated Data – No Correction Factors Applied)

Sample: Indium




Indium melt from 155 °C to 157 °C



Verifying Cell Constant & Temperature

- Run Indium as a sample (i.e. in standard mode not calibration mode)
- Analyze melt and record melt onset & heat of fusion
- Compare to known values
 - Melting of In **156.6°C**
 - Heat of Fusion **28.71J/g**

Verifying Cell Constant & Temperature



Reference Standards for Calibration

- Enthalpy (cell constant)

- Benzoic acid (147.3 J/g) $T_m = 123^{\circ}\text{C}$
- Urea (241.8 J/g) $T_m = 133^{\circ}\text{C}$
- **Indium (28.71 J/g)**
 $T_m = 156.6^{\circ}\text{C}$
- Anthracene (161.9 J/g) $T_m = 216^{\circ}\text{C}$

* GEFTA recommended
Thermochim. Acta, 219 (1993) 333.

ITS 90 Fixed Point

E Zone refined organic compound
(sublimes)

- Temperature

- Cyclopentane* -150.77°C
- Cyclopentane* -135.09°C
- Cyclopentane* -93.43°C
- Cyclohexane # -83°C
- Water # 0°C
- Gallium # 29.76°C
- Phenyl Ether # 30°C
- p-Nitrotoluene E 51.45°C
- Naphthalene E 80.25°C
- **Indium # 156.60°C**
- Tin # 231.95°C
- Lead* 327.46°C
- Zinc # 419.53°C

Traceable Calibration Materials

- NIST DSC calibration materials:
 - SRM 2232 Indium, $T_m = 156.5985^{\circ}\text{C}$
 - SRM 2220 Tin, $T_m = 231.95^{\circ}\text{C}$
 - SRM 2222 Biphenyl, $T_m = 69.41^{\circ}\text{C}$
 - SRM 2225 Mercury, $T_m = -38.70^{\circ}\text{C}$
- NIST: Gaithersburg, MD 20899-0001
 - Phone: 301-975-6776
 - Fax: 301-948-3730
 - Email: SRMINFO@nist.gov
 - Website: <http://ts.nist.gov/srm>

Traceable Calibration Materials

- LGC DSC Calibration Materials:
 - LGC2601: Indium (TA p/n: 915060-901)
 - LGC2608: Lead
 - LGC2609: Tin
 - LGC2611: Zinc
- Laboratory of the Government Chemist, UK
 - Phone: 44 (0) 181 943 7565
 - Fax: 44 (0) 181 943 7554
 - Email: orm@lgc.co.uk


Verifying Baseline

- Run Empty cell (no pans), -90°C to 400°C (w/ RCS) at 20°C/min
 - Experiment is run in the standard mode
 - Plot mW vs. temperature on a 1mW scale
 - Should look fairly flat on this scale
 - Should be around zero heat flow
 - Measure drift and compare to instrument specifications.
 - Verify performance periodically


Verifying Baseline

- Importance of a flat baseline:
 - Detecting very weak transitions
 - Accurate integration of enthalpy
 - Kinetics, partial area analysis, extent of reactions
 - Initial crystallinity


Verifying Baseline



Verifying Baseline



Empty Cell Baseline at 20 °C/min – DSC2500



Experimental Design: Instrument Set Up



Instrument Hardware and Gas Selection Considerations



Temperature Range Dependent On The Cooling System

- Finned Air Cooling System (FACS): Ambient to 725°C
- Quench Cooling Accessory (QCA): -180°C to 400°C
- Liquid Nitrogen Cooling System (LN2P): -180°C to 550°C
- RCS120: -120°C to 400°C
- RCS90: -90°C to 550°C
- RCS40: -40°C to 400°C



Purge Gas Selection

- **Nitrogen**
 - inert, inexpensive and readily available
 - flow rate of 50ml/min
- **Helium**
 - a high thermal conductivity gas which improves response time and cooling capabilities
 - the recommended purge gas when using the LN2 accessory at temperatures below -100°C
 - flow rates of 10-25ml/min are typically used; cell constant affected by flow rate
- **Air/Oxygen**
 - used when studying oxidative stability of materials




Sample Press and Pan Selection

- Aluminum: max. temperature of 600°C
- Gold
- Copper
- Graphite, Alumina
- Platinum
- Stainless Steel

Cooling Accessories

- Finned Air Cooling System (FACS): Ambient to 725°C
- Quench Cooling Accessory (QCA): -180°C to 550°C *
- Liquid Nitrogen Cooling System (LNCS): -180°C to 550°C
- RCS120: -120 °C to 400 °C
- RCS90: -90°C to 550°C
- RCS40: -40°C to 400°C



Performance of the Cooling Accessories

RCS90 Controlled Cooling Rates, from 550°C (upper limit)*

Controlled Rate	To Lower Temperature
100°C/min	300°C
50°C/min	120°C
20°C/min	-20°C
10°C/min	-50°C
5°C/min	-75°C
2°C/min	-90°C

RCS40 Controlled Cooling Rates, from 400°C (upper limit)*

Controlled Rate	To Lower Temperature
65°C/min	250°C
50°C/min	175°C
20°C/min	40°C
10°C/min	0°C
5°C/min	-15°C
2°C/min	-40°C


LNCS Controlled Cooling Rates, from 550°C (upper limit)*

Controlled Rate	To Lower Temperature
100°C/min	200°C
50°C/min	0°C
20°C/min	-100°C
10°C/min	-150°C
5°C/min	-165°C
2°C/min	-180°C

* Performance may vary slightly, depending on laboratory conditions.

Effect of Cooling Rate on Heat Flow


High Density Polyethylene 10C/min Cooling Rate with RCS90



Selecting the Cooler – Discovery DSC

TA Instruments TRIOS

Application
Discovery DSC
Information
General
Cooler
Auto Sampler
Temperature Cal
Heat Capacity

 Cooler Settings

This is used to select the cooler type

Cooler Selection

Cooler Selection RCS 90 Cooler ▼

☒ Activate secondary purge when lid is opened (R)


Between Runs

☒ Leave Cooler On


☒ Autofill LN2P if below %

Finned Cooler
Quench Cooler
LN2P Cooler
RCS 40 Cooler
RCS 90 Cooler

Selecting the Cooler – Q-series



Selecting the Purge Gas – Discovery DSC



Setting the Purge Gas Flow Rate – Discovery DSC

This is used to select which gas is going to the DSC cell and the flow rate for that gas.

General

Temperature: 40.01 °C

Flange Temperature -82.61 °C

Cooler Selection RCS 90 Cooler

Standby 40.0 °C

RCS ON

Gas 1 Nitrogen: 50.0 mL/min

Base purge 263.4 mL/min

Gas 1: Nitrogen

Flow rate 50

Selecting the Purge Gas – Q-series DSC

The screenshot displays the QSeries software interface for a Q1000 DSC. The 'Tools' menu is open, and 'Instrument Preferences...' is selected. The 'Instrument Preferences' dialog box is shown with the 'MFC Purge' tab active. Both 'Gas #1' and 'Gas #2' are set to 'Nitrogen'. The checkbox 'Stop experiment when flow rate deviates from the set value' is checked. A red box highlights the 'Tools' menu and the 'Instrument Preferences...' option. Another red box highlights the 'Instrument Preferences' dialog box. A third red box highlights the text 'This is used to specify the type of gas connected to Gas #1 and Gas #2 inlets'.

This is used to specify the type of gas connected to Gas #1 and Gas #2 inlets

The background shows a graph of Heat Flow (mW) vs. Temperature (°C) with a cooling curve. The y-axis ranges from 0.80 to 2.00 mW, and the x-axis ranges from 0.80 to 2.00 °C. The graph shows a baseline at approximately 1.40 mW, a sharp endothermic peak at approximately 1.60 mW, and a return to the baseline at approximately 1.80 mW. The peak is labeled 'COOLING GAS 1/20 PSIG INLET MAX'.

Setting the Purge Gas Flow Rate – Q-series DSC

Run 1: Standby Temp: 27.50 °C Store: Off Gas: 1 Event: On

Sequence
Run 1

Notes
Operator: LEW
Pan Type: Aluminum
Extended Text:

Mass Flow Control Settings
Sample: #1 - Nitrogen Flow Rate: 50 mL/min

Signal

Signal	Value
Method Time	0.00 min
Segment Time	0.00 min
Remaining Run Time	0 min
Temperature	27.50 °C
Heat Flow	0.037 mW
Heat Capacity	0.000 mJ/°C
Sample Purge Flow	50.01 mL/min
Set Point Temp	27.50 °C
Heater Power	38.568 W
Flange Temperature	-91.59 °C
Heater Temperature	28.15 °C

Running Segment Description

- 1 Equilibrate at -90.00 °C
- 2 Ramp 10.000 °C/min to 300.00 °C

Heat Flow (mW)

Temperature (°C)

01 28.00 min. Append Apply Cancel Help

Ready Stand by Standard Seg 0 in Run 1 13:26:45

This is used to select which gas is going to the DSC cell and the flow rate for that gas.

Recommended Purge Gas Flow Rates

Module

All TA DSC's

Purge Port

50 ml/min (N₂) or 25 ml/min (He)

- If purge gas is too slow - possible moisture accumulation & early aging of the cell
- If purge gas is too fast – excessive noise

Experimental design: Sample Preparation and Considerations



TGA for DSC Experimental Design


- Thermogravimetric Analysis (TGA) measures weight loss or gain as a function of temperature, time and atmosphere.
- General applications of TGA include:
 - thermal stability
 - residual solvent, out gassing, moisture sorption/desorption
 - filler/fiber content
 - weight loss on cure
- TGA measurements are extremely useful in selecting experimental conditions for DSC experiments and for interpreting results.



Selecting Optimum Experimental Conditions

- If possible, run a TGA experiment before beginning DSC tests on new materials
- Heat approximately 10mg sample in the TGA at 10°C/min to determine:
 - Volatile content
 - Unbound water or solvent is usually lost over a broader temperature range at a lower temperature than a hydrate/solvate
 - Decomposition temperature
 - DSC results are of little value once the sample has lost 5% weight due to decomposition (not desolvation)
 - Decomposition is a kinetic process (time & temperature dependent). The measured decomposition temperature will shift to lower temperatures at slower heat rates

Typical TGA data: TGA of Drug A Monohydrate



Selecting Optimum Experimental Conditions

- Use TGA data to help select DSC experimental conditions
 - Standard (non-hermetic) vs. Hermetic (sealed) pans
 - Use hermetic pan if sample loses approximately 0.5% weight or more.
 - Use hermetic pan with pin hole lids if sample loses volatiles such as water
 - Maximum Temperature
 - Excessive decomposition will contaminate the DSC cell between runs
 - When comparing samples, always use the same experimental conditions

DSC Pan Selection

- High thermal conductivity – aluminum, gold, copper, platinum
- Inert - alodined aluminum, ceramic, graphite
- Flatness of the pan for optimal thermal contact
- Standard, non-hermetic vs. hermetic sealing
- Capacity/sample volume
- Temperature range



Sample Pans

- Type of pan depends on:
 - Sample form
 - Volatilization
 - Temperature range
- Use lightest, flattest pan possible
- Always use reference pan of the same type as sample pan

Tzero Press (P/N 901600.901)




Tzero Press kit includes die sets for:

1. Tzero Pans / Tzero Lids and Tzero Low-Mass Pans / Tzero Lids (Black)
2. Tzero Pans / Tzero Hermetic Lids (Blue)
3. Standard Aluminum Pans / Lids (Green)
4. Standard Hermetic Pans / Lids (White)

The kit also includes one box each of Tzero Pans (100) and Tzero Lids (100).

TA Instruments Tzero Pans (Aluminum)

Tzero Pan



- The Tzero pan has been engineered to have a perfectly flat bottom and not to deform during crimping. This ensures the optimal contact between pan and sensor, minimizing the contact resistance and improving resolution.
- The Tzero Pan can be configured for non-hermetic or hermetic use. P/N 901683.901 Tzero Pans (pkg. of 100)


Tzero Low-Mass Pan



- The Tzero Low-Mass Pan is designed for the highest sensitivity when sample mass is limited. P/N 901670.901 Tzero Low-Mass Pans (pkg. of 100). Can only be used with the non-hermetic Tzero lid.


TA Instruments Tzero Pans

Tzero Lid



- Tzero Lid (P/N: 901671.901) - Lightweight aluminum lids for use in sample encapsulation with the Tzero Pans and the Tzero Low-Mass Pans. The seal is not hermetic.


Tzero Hermetic Lid



- Tzero Hermetic Lid (P/N: 901684.901) (pkg. of 100) and P/N: 901685.901 Tzero Hermetic Pinhole Lid (75 micron diameter pinhole) (pkg. of 50). Used only with the Tzero pan, not the low mass Tzero pan

Standard Series DSC Aluminum Pans

- Part numbers for the pans and lid
 - 900760.901 Classic Aluminum Pans (pkg. of 200) (higher sidewall compared to the standard aluminum sample pans to accommodate larger samples)
 - 900786.901 Aluminum Sample Pans (pkg. of 200)
 - 900779.901 Aluminum Lids (pkg. of 200)
- Pan & lid weighs ~23mg, bottom of pan is flat
- use up to **600°C maximum**
- Used for solid non-volatile samples
- Always use lid (see exceptions)
 - Lid improves thermal contact
 - Keeps sample from moving
- Exceptions to using a lid
 - Running oxidative experiment
 - Running PCA experiment



Standard Series DSC Other Pans (Non Hermetic)


- Standard pans are available in:
 - Gold (p/n 900866.901 pan, p/n 900868.901 lid): up to 725°C
 - Graphite (p/n 900874.901 pan, p/n 900873.901 lid): up to 725°C (in N₂)
- Standard Pans with no lids available
 - Platinum (p/n 900578.901): up to 725°C
 - Copper (p/n 900867.901) : up to 725°C

Sample Shape

- Keep sample thin
- Cover as much the bottom of pan as possible



Hermetic DSC Pans



Aluminum
Hermetic
3 atm



Sample:
1. Liquid
2. Solid with
volatile
content

High Volume
Stainless
Steel
40 atm





Gold
Hermetic
6 atm





Standard Series Aluminum Hermetic Pans

- Part numbers for the pans and lid
 - 900793.901 Aluminum Sample Pans, Hermetic (pkg. of 200)
 - 900794.901 Aluminum Lids, Hermetic (pkg. of 200)
 - 900860.901 Hermetic Lids with Pinhole (pkg. of 50)
- Pan & Lid weigh ~55mg, bottom of pan is not as flat as standard pans
- Used for liquid samples and samples with volatiles
- Always use lid (same exceptions as before)
- After sealing pans, the lid should form a dome




Hermetic Pans Summary

- Hermetic Pans are available in:
 - Aluminum: $<600^{\circ}\text{C}$; <3 atm (300 kPa gauge)
 - Alodined Aluminum: $<200^{\circ}\text{C}$; <3 atm (300 kPa gauge)
 - Gold: $<725^{\circ}\text{C}$; <6 atm (600 kPa gauge)
- Specialized Sealed Pans
 - High Volume: $100\mu\text{L}$; $<250^{\circ}\text{C}$; 37 atm (575 psi)
 - P/N 900825.901
 - High Pressure: $35\mu\text{L}$; $<300^{\circ}\text{C}$; 100 atm (1450 psi)
 - P/N 900808.901



Note: 3 atm is approximately 44 psi

It Can Matter What Pan You Use



What if Sample Spills out of the Pan?

Keeping the DSC Cell Clean

- One of the first steps to ensuring good data is to keep the DSC cell clean
- How do DSC cells get dirty?
 - Decomposing samples during DSC runs
 - Samples spilling out of the pan
 - Transfer from bottom of pan to sensor

Cleaning the Cell

- Use solvent – slightly damp swab with an appropriate solvent
 - Heat cell to 200°C for 10 min to drive off any remaining solvent
- If the cell is still dirty
 - Clean w/ brush
 - Be careful with the Tzero™ thermocouple
 - Fibers in cell from cleaning brush need to be removed



Cleaning Cell: Bake Out Procedure

- Bake out
 - Should be used as a last resort if none of the previous steps are effective
 - Involves Air purge and/or an open lid
 - Heat at 20°C/min to appropriate temp (max of 550°C on Q series, max. 400°C in Discovery)
 - Do NOT hold Isothermal at the upper temperature
 - Cool back to room temp & brush cell again
- Irrespective of the cleaning method used, always verify the baseline at the end of the cleaning procedure, and recalibrate the DSC if required
- Check out the TA Tech tip video on cleaning the DSC cell:
<https://www.youtube.com/watch?v=cclJXrbUICA>

Sample Preparation: Optimization of Sample Mass

- Sample Preparation
 - Weight of 5-10 mg for polymers; 10-15 mg for cross-linked thermosets; 3-5 mg for metal or chemical melting
 - Goal is to achieve a change of 0.1-10mW heat flow in going through the transition

$$\frac{dH}{dt} = C_p \frac{dT}{dt} + f(T, t)$$

Heat Flow Change During a Transition

